



PRODUCT  
CATALOGUE

2019-2020

# TABLE OF CONTENTS

---

Introduction .....	1
Vision, Objectives, Mission .....	1
Sites, Our Products .....	2
Standards, Our Markets .....	3
MILESTONES .....	4
LATEX KITS .....	5
TURBIDIMETRIC LATEX KITS .....	6
SYPHILIS KITS .....	6
FEBRILE ANTIGENS .....	7
BLOOD GROUPING REAGENTS.....	8
HEMATOLOGY TESTS .....	9
SICKLE CELL KITS .....	9
INFECTIOUS DISEASE RAPID TESTS .....	10
ANTIBODY TESTING .....	10
ANTIGEN TESTING .....	11
URINE REAGENT STRIPS .....	12
FERTILITY RAPID TESTS .....	13
KIDNEY FUNCTION RAPID TESTS .....	13
INFLAMMATION AND CANCER MARKERS .....	14
CARDIAC MARKERS RAPID TESTS .....	14
DOA RAPID TESTS .....	15
CLINICAL CHEMISTRY KITS .....	17
ALCOHOL TESTS .....	18
STAINS FOR HISTOLOGY & MICROBIOLOGY .....	19
ANTIBIOTIC SENSITIVITY DISCS .....	20
ELISA KITS .....	21
HORMONE ELISA KITS .....	21
TORCH ELISA KITS .....	22
INFECTIOUS DISEASES ELISA KITS .....	22
OTHER ELISA KITS .....	22
HOME TESTS .....	23
BLOOD GLUCOSE MONITORING SYSTEMS .....	24
CERTIFICATES .....	25
CONTACT US .....	26
INTERNATIONAL PRESENCE .....	27

# INTRODUCTION

Atlas Medical was established in 1996 as a manufacturer and supplier of quality Diagnostic Reagents and Kits. Our products are sold in over 80 countries worldwide.

The company is located at the Cambridge Science Park, Cambridge, UK. In addition to the UK site, the company has offices in Germany and Turkey as well as two purpose-built modern facilities in both Jordan and Malaysia. We take quality assurance very seriously and strive to produce goods to the highest standards known in the industry, including, ISO13485 & CE mark and US FDA standards. Our R&D team constantly develops and innovates novel products that significantly contribute to the advancement of the Diagnostic Industry.

## Vision

To be a major provider of quality medical diagnostic products to local, regional and international markets.

## Mission

Our mission is to develop, produce and provide our customers with high quality products and excellent customer services through deep understanding of customers' needs and perception, recruitment of high caliber professionals & technicians, adopting strict quality assurance and control procedures and embracing new scientific advancements in the medical lab diagnostic field.

## Objectives



*High and Consistent Quality*



*Satisfied Customer*



*Continuous Improvement & Innovation*



*Care for the Environment & Working Conditions*



Ludwig-Erhard-Ring 3  
15827 Blankenfelde-Mahlow,  
Berlin, Germany  
Tel: +49 - 33708 - 3550 30  
Email: amug@atlas-medical.com



William James House, Cowley  
Road, Cambridge, CB4 0WX, UK.  
Tel.: +44 (0) 1223 858 910  
Fax: +44 (0) 1223 858 524  
Email: info@atlas-medical.com



Unit G-E-2A, Enterprise 4,  
Technology Park Malaysia,  
Lebuhraya Puchong -  
Sungai Besi, 57000  
Bukit Jalil, Kuala Lumpur.



Sahab Free Zone Area,  
Amman 11512 Jordan  
Tel: +962 6 4026468  
Fax +962 6 4022588  
Email : info@atlas-medical.com



Yenigün Mahallesi Bağcılar  
Cad. ünlü Apt. Apt. NO:131/1  
Bağcılar / Istanbul / Turkey  
Tel: +90 537 796 92 63  
Email:info@atlas-medical.com

## Our Products

---

- 🌀 Latex related kits, indirect hemo-agglutination, carbon antigens (serology)
- 🌀 Blood grouping and rare sera (blood banking)
- 🌀 Coagulation and D-Dimer reagents (hematology)
- 🌀 Stained bacterial antigens, blood culture bottles (febrile antigens - microbiology).
- 🌀 Rapid tests and urine reagent strips (lateral flow immunoassay and dry chemistry).
- 🌀 Stains for microbiology and histopathology
- 🌀 ELISA kits (infectious diseases, hormones, torch and Vitamins)
- 🌀 Clinical chemistry (biochemistry)
- 🌀 Drug-of-Abuse and alcohol tests (Toxicology)
- 🌀 Home (OTC) Tests: Fertility Home Tests & General Health Screening Tests.
- 🌀 Blood Glucose Monitoring Systems
- 🌀 Quantas Lateral Flow



## Standards



Our products are manufactured in accordance to the standards as set in the European In-Vitro Diagnostic Directive 98/79/EC. This has led to the successful attainment of Annex IV Full Quality Assurance Certification and the declaration of conformity for CE marking purposes for many of our IVD products, either self-declared or through our Notified Body LNE/G-MED in France.



To complete the quality assurance scheme the company has put in place a robust Quality Management and Enhancement System that has concluded in the successful attainment of ISO13485: 2016 certificate




the company also adheres to the US-FDA regulations and had already FDA-cleared few products for the US market. Our products are registered in numerous countries.



## Our Markets

Atlas Medical enjoys a good presence in many international markets. We take pride in our export activities through our dedicated export department. We actively participate in major industry-related exhibitions seeking keen representatives around the globe to sell and distribute our products in their respective countries. We are internationally represented in more than 80 of countries spanning in five continents: Europe, North America, South America, Africa and Asia. Our efforts will continue to increase our representation to include most markets around the globe.

# MILESTONES



# LATEX KITS



## Overview

Latex kits offer a quick and simple assay to diagnose a range of pathogens and medical conditions. The assay is based on an immunological reaction between the detected analyte in the sample and its corresponding antibody or antigen already coated on latex particles.

## Features

- 🔗 They cover a selection of routine tests in serology and microbiology.
- 🔗 They are conveniently packed in sizes of 50 or 100 tests and includes all the necessary reagents, controls, stirrers and slides to conduct the test.
- 🔗 Affordable, easy to use, dependable and offer a clear and visible agglutination for doubt-free results within 2 minutes.
- 🔗 Some Latex Kits come with a Buffer.

Item Code	Item Description	Sizes
8.00.00.0.0050 8.00.00.0.0100	CRP Latex Kit	50 Tests 100 Tests
8.00.01.0.0050 8.00.01.0.0100	CRP Latex Kit with Buffer	50 Tests 100 Tests
8.00.02.0.0050 8.00.02.0.0100	ASO Latex Kit	50 Tests 100 Tests
8.00.03.0.0050 8.00.03.0.0100	ASO Latex Kit with Buffer	50 Tests 100 Tests
8.00.04.0.0050 8.00.04.0.0100	RF Latex Kit	50 Tests 100 Tests
8.00.05.0.0050 8.00.05.0.0100	RF Latex Kit with Buffer	50 Tests 100 Tests
8.00.07.0.0050 8.00.07.0.0100	hCG Latex Kit	50 Tests 100 Tests
8.00.11.0.0050 8.00.11.0.0100	SLE Latex Kit	50 Tests 100 Tests
8.00.16.0.0050 8.00.16.0.0100	Rota Virus Latex Kit	50 Tests 100 Tests
8.00.17.0.0050 8.00.17.0.0100	D-Dimer Latex Kit	50 Tests 100 Tests
8.00.21.0.0050 8.00.21.0.0100	Waler Rose Kit	50 Tests 100 Tests
8.00.31.0.0050 8.00.31.0.0100	IM Latex Kit	50 Tests 100 Tests



Item Code	Item Description	Sizes
8.00.09.0.0050 8.00.09.0.0100	Toxo Latex Kit	50 Tests 100 Tests
8.00.10.0.0050 8.00.10.0.0100	Toxo Latex Kit with Buffer	50 Tests 100 Tests
8.00.14.0.0100	Rubella Latex Kit	50 Tests 100 Tests

# TURBIDIMETRIC LATEX KITS

## Overview

The turbidimetric assay is based on the agglutination reaction between latex particles coated with antibody and the antigen in solution. The intended use for Turbilatex products is to detect and quantify the antigen present in human serum samples.



## Features

- Atlas Medical offers a dynamic range of Turbidimetric Latex Kits which are conveniently packed in sizes of 50, 100 and 250 tests and include all the necessary accessories.

Item Code	Item Description	Sizes
8.44.00.0.0050	RF Turbidimetric Latex Kit	50 Tests
8.44.00.0.0100		100 Tests
8.44.00.0.0250		250 Tests
8.44.01.0.0050	CRP Turbidimetric Latex Kit	50 Tests
8.44.01.0.0100		100 Tests
8.44.01.0.0250		250 Tests
8.44.02.0.0050	ASO Turbidimetric Latex Kit	50 Tests
8.44.02.0.0100		100 Tests
8.44.02.0.0100		250 Tests
8.44.03.0.0050	D-Dimer Turbidimetric Latex Kit	50 Tests
8.44.03.0.0100		100 Tests
8.44.03.0.0250		250 Tests

# SYPHILIS KITS

## Overview

Atlas Medical offers a number of assays to detect Syphilis that include: TPHA kits which are used for the detection of antibodies to *Treponema pallidum* in human Serum using micro haemagglutination; VDRL and RPR kits which are based on non-Treponemal flocculation to detect reagin antibodies in serum or plasma; and rapid lateral flow immunoassay tests to detect Syphilis antibodies in serum, plasma or whole blood.



Item Code	Item Description	Sizes
8.00.18.0.0100	RPR Carbon Antigen Kit	100 Tests
8.00.18.0.0250		250 Tests
8.00.18.0.0500		500 Tests
8.00.19.0.0050	TPHA Kit	50 Tests
8.00.19.0.0100		100 Tests
8.00.19.0.0200		200 Tests
8.00.20.0.0250	VDRL Kit	250 Tests
8.00.20.0.0500		500 Tests
8.00.20.0.2500		2500 Tests

## Features

- Easy to use, affordable and conveniently packed in different sizes to suit all needs.
- They include all the necessary reagents/devices, controls, stirrers and slides to conduct the test





# FEBRILE ANTIGENS



## Overview

Febrile antigen kits are based on bacterial suspensions that agglutinate in the presence of antibodies formed in human infection by certain fever-causing microbial agents. In positive samples, the agglutination is macroscopically visible at certain antibody levels in serum. These antigen reagents are used for the qualitative and semi quantitative febrile screening purposes.

## Features

- Atlas Medical Febrile Antigen kits contain various types of antigens for Brucella, Proteus, Salmonella typhi and paratyphi, and their controls as needed.
- Atlas Medical Febrile Antigen kits are competitively priced and easy to use, and give clear results within 2 minutes



Item Code	Item Description	Sizes
8.01.00.0.0005 8.01.00.0.0050 8.01.00.0.0100	Brucella Rose Bengal Kit	5ml/vial 50 Tests 100 Tests
8.01.01.0.0005 8.01.01.1.0040	Salmonella OA Reagent	5 ml/vial 5x8 ml
8.01.02.0.0005 8.01.02.1.0040	Salmonella OB Reagent	5 ml/vial 5x8 ml
8.01.03.0.0005 8.01.03.1.0040	Salmonella OC Reagent	5 ml/vial 5x8 ml
8.01.04.0.0005 8.01.04.1.0040	Salmonella OD Reagent	5 ml/vial 5x8 ml
8.01.05.0.0005 8.01.05.1.0040	Salmonella HA Reagent	5 ml/vial 5x8 ml
8.01.06.0.0005 8.01.06.1.0040	Salmonella HB Reagent	5 ml/vial 5x8 ml
8.01.07.0.0005 8.01.07.1.0040	Salmonella HC Reagent	5 ml/vial 5x8 ml
8.01.08.0.0005 8.01.08.1.0040	Salmonella HD Reagent	5 ml/vial 5x8 ml
8.01.10.0.0005 8.01.10.1.0040	Brucella Abortus Reagent	5 ml/vial 5x8 ml
8.01.11.0.0005 8.01.11.1.0040	Brucella Melitensis Reagent	5 ml/vial 5x8 ml
8.01.12.0.0005 8.01.12.1.0040	Proteus OX2 Reagent	5 ml/vial 5x8 ml
8.01.13.0.0005 8.01.13.1.0040	Proteus OX19 Reagent	5 ml/vial 5x8 ml
8.01.14.0.0005 8.01.14.1.0040	Proteus OXK Reagent	5 ml/vial 5x8 ml
8.01.15.0.0010	Brucella Antigen Kit ( Brucella melitensis, Brucella abortus )	5 ml/vial
8.01.15.1.0010	Brucella Antigen Kit with Controls (Brucella melitensis, Brucella abortus, 2x1.0 ml Controls )	5 ml/vial
8.01.16.0.0040	Salmonella Antigen Set ( 8 Antigens: OA, OB, OC, OD, HA, HB, HC, HD )	5x8 ml
8.01.16.1.0040	Salmonella Antigen Set ( 8 Antigens: OA, OB, OC, OD, HA, HB, HC, HD ) with 2x1.0 ml Controls	5x8 ml
8.01.17.1.0050	Febrile Antigen Set (10 Antigens: Salmonella OA, OB, OC, OD, HA, HB, HC, HD, Brucella abortus, melitensis) with 2x1.0ml Controls	5x10 ml
8.01.18.0.0030	Salmonella Antigen Set, Widal Kit (6 Antigens: OA, OB, OD, HA, HB, HD)	5x6 ml
8.01.18.1.0030	Salmonella Antigen Set, Widal Kit (6 Antigens: OA, OB, OD, HA, HB, HD) with 2x1.0ml Controls	5x6 ml
8.01.14.0.0005 8.01.14.1.0040	Febrile Antigens Positive Control	1 ml/vial 5 ml/vial
8.01.14.0.0005 8.01.14.1.0040	Febrile Antigen Negative Control	1 ml/vial 5 ml/vial

# BLOOD GROUPING REAGENTS

## Overview

Blood Grouping reagents are used for the identification of blood types. The test procedure is based on the agglutination principle, where red cells possessing the typing antigen agglutinate in the presence of the corresponding antibody in the testing reagent indicating the presence of the tested antigen. No agglutination indicates the absence of the tested antigen.



## Features

- Atlas Medical ABO reagents are prepared from In-Vitro culture supernatants of hybridized immunoglobulin-secreting mouse cell lines.
- The reagents are formulated and optimized for use in tube and slide methods.
- Atlas Medical provides high quality blood grouping reagents that are accurate, easy to use, competitively priced, and conveniently packed in different sizes and options.



Item Code	Item Description	Sizes
8.02.14.0.0010	Anti-D Monoclonal (IgM), Clone 1, 10ml/vial	10 ml/vial
8.02.15.0.0010	Anti-D Monoclonal (IgM), Clone 2, 10ml/vial	10 ml/vial
8.02.16.0.0005	Anti-A1, Lectin (Dolichosbiflorus), 5ml/vial	5 ml/vial
8.02.17.0.0002	Anti-H, Lectin (Ulexeuropaeus), 2ml/vial	2 ml/vial
8.02.18.0.0005	Anti-C Monoclonal, 5ml/vial	5 ml/vial
8.02.19.0.0005	Anti-c Monoclonal, 5ml/vial	5 ml/vial
8.02.20.0.0005	Anti-E Monoclonal, 5ml/vial	5 ml/vial
8.02.21.0.0005	Anti-e Monoclonal, 5ml/vial	5 ml/vial
8.02.22.0.0005	Anti-C+D+E Monoclonal, 5ml/vial	5 ml/vial
8.02.27.0.0002	Anti-Fya, Human, 2ml/vial	2 ml/vial
8.02.28.0.0002	Anti-Fyb, Human, 2ml/vial	2 ml/vial
8.02.29.0.0002	Anti-k, Human, 2ml/vial	2 ml/vial
8.02.30.0.0002	Anti-Kpa, Human, 2ml/vial	2 ml/vial
8.02.31.0.0002	Anti-Kpb, Human, 2ml/vial	2 ml/vial
8.02.32.0.0002	Anti-Jka, Human, 2ml/vial	2 ml/vial
8.02.36.0.0005	Anti-K Monoclonal, 5ml/vial	5 ml/vial



Item Code	Item Description	Sizes
8.02.00.0.0010 8.02.00.1.0100	Anti-A Monoclonal reagent (titer: 1/512)	10 ml/vial 10x10 ml
8.02.01.0.0010 8.02.01.1.0100	Anti-B Monoclonal Reagent (Titer: 1/512)	10 ml/vial 10x10 ml
8.02.02.0.0010 8.02.02.1.0100	Anti-AB Monoclonal Reagent (Titer: 1/512)	10 ml/vial 10x10 ml
8.02.03.0.0010 8.02.03.1.0100	Anti-D IgG/IgM Blend Reagent (Titer: 1/128)	10 ml/vial 10x10 ml
8.02.04.0.0010 8.02.04.0.0100	Anti-A Monoclonal Reagent (Titer: 1/256)	10 ml/vial 10x10 ml
8.02.05.0.0010 8.02.05.0.0100	Anti-B Monoclonal Reagent (Titer: 1/256)	10 ml/vial 10x10 ml
8.02.06.0.0010 8.02.06.1.0100	Anti-AB Monoclonal Reagent (Titer: 1/256)	10 ml/vial 10x10 ml
8.02.07.0.0010 8.02.07.1.0100	Anti-D IgG/IgM Blend Reagent (Titer: 1/64)	10 ml/vial 10x10 ml
8.02.08.0.0010 8.02.08.1.0100	Bovine Albumin 22%	10 ml/vial 10x10 ml
8.02.09.0.0010 8.02.09.1.0100	Bovine Albumin 30%	10 ml/vial 10x10 ml
8.02.10.0.0010 8.02.10.1.0100	Anti-Human Globulin (Green) (Titer 1/512)	10 ml/vial 10x10 ml
8.02.11.0.0010 8.02.11.1.0100	Anti-Human Globulin (Green) (Titer 1/256)	10 ml/vial 10x10 ml
8.02.47.0.0030	ABO Set (Anti-A (1/512), Anti-B (1/512), Anti-D (1/128))	3x10 ml
8.02.47.1.0030	ABO Set (Anti-A (1/265), Anti-B (1/265), Anti-D (1/64))	3x10 ml
8.02.49.0.0040	ABO Set (Anti-A (1/256), Anti-B (1/256), Anti-AB (1/256), Anti-D (1/64))	4x10 ml
8.02.53.0.0040	ABO Set (Anti-A (1/512), Anti-B (1/512), Anti-AB (1/512) Anti-D (1/128))	4x10 ml
8.02.52.0.0010 8.02.52.0.0100	Rh-D Negative Control	10 ml/vial 10x10 ml
8.02.63.1.0010 8.02.63.1.0100	Antibody Enhancement Solution (LISS)	10 ml/vial 10x10 ml
8.02.23.0.0002	Anti-M, Human, 2ml/vial	2 ml/vial
8.02.24.0.0002	Anti-N, Lectin (Viciagraminea), 2ml/vial	2 ml/vial
8.02.25.0.0002	Anti-S, Human, 2ml/vial	2 ml/vial
8.02.26.0.0002	Anti-s, Human, 2ml/vial	2 ml/vial
8.02.34.0.0002	Anti-Lua, Human, 2ml/vial	2 ml/vial
8.02.37.0.0002	Anti-Lea, Monoclonal, 2ml/vial	2 ml/vial
8.02.38.0.0002	Anti-Leb, Monoclonal, 2ml/vial	2 ml/vial
8.02.39.0.0002	Anti-P1, Monoclonal, 2ml/vial	2 ml/vial

# HEMATOLOGY TESTS

## Overview

Atlas Medical supplies coagulation reagents and D-Dimer tests. The coagulation reagents include PT and PTT in liquid formats and in various sizes to suit most lab applications. The range also includes normal and abnormal coagulation controls. Atlas Medical D-Dimer test is based on latex agglutination and is perceived as one of the highest quality.

## Features

- The kit includes positive & negative controls and assay accessories.
- The Kit comes in sizes of 50 and 100 tests.
- Atlas Medical provides high quality blood grouping reagents that are accurate, easy to use, competitively priced, and conveniently packed in different sizes and options.



### Co-aggulation Reagents

Item Code	Item Description	Sizes
8.02.40.1.0010 8.02.40.1.0050 8.02.40.1.0100	PT Calcium Rabbit Brain Thromboplastin, Liquid	2ml (20 Tests) 5ml (50 Tests) 10ml (100 Tests)
8.02.41.1.0040 8.02.41.1.0100 8.02.41.1.0200	APTT (PTT) Micronised Silica Platelet Substitute	2ml (40 Tests) 5ml (100 Tests) 10ml (200 Tests)
8.02.44.0.0040 8.02.44.0.0100	PT Kit with Normal Control	2x2ml + 1ml 2x10ml + 1ml
8.02.45.0.0080 8.02.44.0.0100	APTT ( PTT) Kit with Normal Control	2x2ml + 1ml 2x10ml + 1ml
8.02.48.0.0010	Calcium Chloride, 25 mM	10ml/vial
8.02.60.0.0006	Normal Coagulation Control	6x1ml
8.02.61.0.0006	Abnormal Coagulation Control	6x1ml
8.02.69.0.0005	Fibrinogen Reagent	5ml/vial



### Hemoglobin Reagents

Item Code	Item Description	Sizes
8.02.46.1.0500 8.02.46.1.1000 8.02.46.1.3000	Drabkins Reagent, 40x ( White Plastic Bottle )	50ml/vial 2x50ml 6x50ml
8.02.44.0.0040 8.02.44.0.0100	PT Kit with Normal Control	2x2ml + 1ml 2x10ml + 1ml



# SICKLE CELL KITS

## Overview

Sickle cell disease (also called sickle cell anemia) is an inherited blood disorder that affects red blood cells. The sickle cell gene causes the body to produce abnormal hemoglobin.

## Features

- Atlas Sickle Cell Kits is a qualitative solubility test for Sickle Haemoglobin.
- The test can be performed in two ways:
  - A screening test to detect sickle haemoglobin (HbS)
  - A centrifugation test to differentiate the sickle cell trait (AS) from sickle cell anaemia (SS).



Item Code	Item Description	Sizes
8.02.67.0.0100	Sickle Cell Kit, 100 Tests	4x50ml Buffer, 4 vials of Sodium Dithionate
8.02.68.0.0001	Sickle Cell positive & negative control set	1ml each

# INFECTIOUS DISEASE RAPID TESTS

## ANTIBODY TESTING

### Overview

Atlas Medical offers an extensive range of lateral flow immunoassay tests for the rapid detection of antibodies and antigens in human samples (blood, serum, plasma, urine, oral swabs, nasal swabs, and feces). This range includes tests to detect a wide variety of viruses, microorganisms and parasites.



### Features

- Atlas Medical infectious disease rapid tests are reliable, accurate and supplied in both cassette and strip formats.
- The kits are conveniently packed in different sizes of 20, 30, 40, 50 and 100 tests per kit and include the necessary test accessories to perform the assay.
- Atlas Medical infectious disease rapid tests can be used with Quantas™ to obtain quantitative results.



### Quantas™

- Quantas™ reader is made with German technology for reliable, fast and accurate results.
- A printer can be attached to document the results.



Item Code	Item Description	Sizes
8.16.96.0.0020	CRP Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	20 Tests/Box
8.46.00.0.0020	Hs-CRP Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	20 Tests/Box
8.16.95.0.0020	D-Dimer Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	20 Tests/Box



Item Code	Item Description	Sizes
8.04.20.0.0001 8.04.20.0.0020	H.pylori Antibody Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.21.0.0001 8.04.21.0.0020	H.pylori Antibody Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.22.0.0001 8.04.22.0.0100	H.pylori Antibody Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box
8.04.41.0.0001 8.04.41.0.0020	Syphilis Antibody Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.42.0.0001 8.04.42.0.0020	Syphilis Antibody Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.43.0.0001 8.04.43.0.0100	Syphilis Antibody Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box
8.16.16.0.0001 8.16.16.0.0020	TB Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box





Item Code	Item Description	Sizes
8.04.27.0.0001 8.04.27.0.0020	HIV 1/2 Antibody Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.28.0.0001 8.04.28.0.0020	HIV 1/2 Antibody Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.29.0.0001 8.04.29.0.0100	HIV 1/2 Antibody Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box
8.04.30.0.0001 8.04.30.0.0020	HCV Antibody Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.31.0.0001 8.04.31.0.0020	HCV Antibody Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.32.0.0001 8.04.32.0.0100	HCV Antibody Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box
8.04.35.0.0001 8.04.35.0.0020	HBs Antibody Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.36.0.0001 8.04.36.0.0100	HBs Antibody Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box

# INFECTIOUS DISEASE RAPID TESTS

## ANTIGEN TESTING



STOOL SAMPLE		
Item Code	Item Description	Sizes
8.04.23.1.0020	H.pylori Antigen Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.04.24.1.0025	H.pylori Antigen Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.04.69.0.0020	Rotavirus Antigen Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.04.70.0.0025	Rotavirus Antigen Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.04.71.0.0020	Adenovirus Antigen Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.04.72.0.0025	Adenovirus Antigen Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.04.73.0.0020	Rota-Adeno Antigens Combo test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.04.74.0.0025	Rota-Adeno Antigens Combo test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.16.01.0.0020	Crypto Virus Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.16.02.0.0025	Crypto Virus Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.16.31.0.0020	Giardia Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.16.30.0.0025	Giardia Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.16.33.0.0020	Crypto-Giardia Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.16.32.0.0025	Crypto-Giardia Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.16.41.0.0020	E.coli Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.16.40.0.0025	E.coli Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.16.82.0.0025	Salmonella typhi Antigen Test Cassette, Stool Sample, Individually Pouched	25 Tests/Box
8.16.85.0.0025	Salmonella paratyphi Antigen Test Cassette, Stool Sample, Individually Pouched	25 Tests/Box
8.16.80.0.0020	Clostridium difficile Antigen Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.16.81.0.0025	Clostridium difficile Antigen Test Strip, Stool Sample, Individually Pouched	25 Tests/Box
8.16.86.0.0025	C.Difficile Toxin A+B, Test Cassette, Stool Sample, Individually Pouched	25 Tests/Box
8.16.91.0.0025	Norovirus Genogroups I & II Ag, Test Cassette, Stool Sample, Individually Pouched	25 Tests/Box



SWAB & NASAL SAMPLE		
Item Code	Item Description	Sizes
8.04.25.0.0020	Strep A Test Cassette, Swab Sample	20 Tests/Box
8.45.00.0.0020	Strep B Test Cassette, Swab Sample	20 Tests/Box
8.45.01.0.0020	Strep A+B Test Cassette, Swab Sample	20 Tests/Box
8.04.86.0.0020	Influenza A+B Test Cassette, Nasal Sample	20 Tests/Box
8.04.96.0.0025	Influenza A+B Test Strip, Nasal Sample	25 Tests/Box
8.16.20.0.0020	RSV Test Cassette, Swab Sample	20 Tests/Box
8.16.22.0.0025	RSV Test Strip, Swab Sample	25 Tests/Box
8.16.37.0.0020	Adeno Respiratory Antigen Test Cassette, Swab Sample	20 Tests/Box
8.16.36.0.0025	Adeno Respiratory Antigen Test Strip, Swab Sample	25 Tests/Box
8.16.39.0.0020	Adeno - RSV Respiratory Test Cassette, Swab Sample	20 Tests/Box
8.16.38.0.0025	Adeno - RSV Respiratory Test Strip, Swab Sample	25 Tests/Box




BLOOD SAMPLE		
Item Code	Item Description	Sizes
8.04.37.0.0020	Malaria Pf. Test Cassette, Whole Blood, Individually Pouched	20 Tests/Box
8.16.14.0.0020	Malaria Pf/Pv. Test Cassette, Whole Blood, Individually Pouched	20 Tests/Box



Item Code	Item Description	Sizes
8.04.26.0.0020	Chlamydia Test Cassette, Urine or Swab	20 Tests/Box
8.16.24.0.0001 8.16.24.0.0020	HBsAg Test Cassette ( Whole Blood/Serum/Plasma), Individually Pouched	Bulk 20 Tests/Box
8.04.33.0.0001 8.04.33.0.0020	HBsAg Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.34.0.0001 8.04.34.0.0100	HBsAg Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box


POSITIVE



NEGATIVE



INVALID



# URINE REAGENT STRIPS

## Overview

Urine Reagent Strips (URS) are widely used in Urinalysis to determine pathological changes in urine. The strips contain dry-chemistry pads that, when dipped in urine, change their colors. The color change allows for the semi-quantitative measurement of various urine parameters. The strips are suitable for lab, point-of-care and even home use.



## Features

- 🔄 Atlas Medical Urine Reagent Strips can be used to detect up to 14 urine parameters.
- 🔄 They are simple to use and the results are visually read within a minute.
- 🔄 The strips are packed in desiccated bottles of 50 or 100 strips.
- 🔄 Atlas Medical can also provide suitable readers to read the strip colors and document the results.



## URINE ANALYZER

Atlas Urine Analyzer is a manual analyzer that detects Photosensitive Diode using the method of Reflectance Photometry. Test Categories include Routine, STAT and QC. Atlas Urine Analyzer has an Automatic calibration for accurate results and easy operation. It can read strips with up to 14 parameters, including Microalbumin/Creatinine/Calcium. It has an option to print results for quick and simple record management.

- 🔄 Accurate.
- 🔄 Reliable.
- 🔄 Convenient.
- 🔄 Easy Data Management.

Item Code	Item Description	Sizes
8.03.00.0.0100	URS 1 Parameter: Glucose	100 Strips
8.03.01.0.0100	URS 1 Parameter: Protein	100 Strips
8.03.02.0.0100	URS 1 Parameter: Ketone	100 Strips
8.03.03.0.0100	URS 2 Parameters: Glucose, Ketone	100 Strips
8.03.04.0.0100	URS 2 Parameters: Glucose, Protein	100 Strips
8.03.06.0.0100	URS 3 Parameters: Protein, pH, Glucose	100 Strips
8.03.07.0.0100	URS 3 Parameters: Glucose, Protein, Ketone	100 Strips
8.03.11.0.0100	URS 4 Parameters: Protein, pH, Specific Gravity, Glucose	100 Strips
8.03.12.0.0100	URS 4 Parameters: Protein, pH, Blood, Glucose	100 Strips
8.03.13.0.0100	URS 5 Parameters: Glucose, Protein, Ketone, pH, Blood	100 Strips
8.03.14.0.0100	URS 6 Parameters: Leukocytes, Nitrite, Protein, pH, Blood, Glucose	100 Strips
8.03.44.0.0100	URS 7 Parameter: Glucose, Ketone, Protein, pH, Blood, Bilirubin, Urobilinogen	100 Strips
8.03.23.0.0100	URS 8 Parameters: Glucose, Protein, pH, Ketone, Urobilinogen, Bilirubin, Blood, Nitrite	100 Strips
8.03.15.0.0100	URS 9 Parameters: Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose	100 Strips
8.03.16.0.0100	URS 10 Parameters: Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose	100 Strips
8.03.18.0.0100	URS 11 Parameters: Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose, Ascorbic Acid	100 Strips
8.001.U120	Urine Analyzer For Clinics, U120	1 Unit
8.002.U500	Urine Analyzer For Hospitals, U500	1 Unit



# FERTILITY RAPID TESTS


## Overview

Atlas Medical Fertility Rapid Tests are based on lateral flow immunoassay for the detection of human chorionic gonadotropin (hCG), Ovulation (LH), and Human Follicular Stimulating Hormone (FSH) in urine. Each of the three tests comes in strip, cassette, or midstream formats and are conveniently packed in sizes to suit lab, point-of-care and home uses.



## Features

- 🔄 Accurate.
- 🔄 Convenient.
- 🔄 Easy to use (add or dip in urine).
- 🔄 Competitively priced.
- 🔄 Results are obtained in 1 to 5 minutes.
- 🔄 Different strip sizes are available.



Item Code	Item Description	Sizes
8.04.00.0.0001 8.04.00.0.0020	hCG Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.01.0.0001 8.04.01.0.0020	hCG Test Cassette, Urine/Serum, Individually Pouched	Bulk 20 Tests/Box
8.04.04.0.0001 8.04.04.0.0100	hCG Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.10.0.0001 8.04.10.0.0100	hCG Test Strip, Urine/Serum, Individually Pouched	Bulk 100 Tests/Box
8.04.13.0.0001 8.04.13.0.0015	hCG Midstream Test, Individually Pouched	Bulk 15 Tests/Box
8.04.14.0.0001 8.04.14.0.0020	LH Test Cassette, Urine, Individually pouched	Bulk 20 Tests/Box
8.04.15.0.0001 8.04.15.0.0100	LH Test Strip, Urine, Individually pouched	Bulk 100 Tests/Box
8.04.16.0.0001 8.04.16.0.0015	LH Midstream Test, Individually Pouched	Bulk 15 Tests/Box
8.04.17.0.0001 8.04.17.0.0020	FSH Test Cassette, Urine, Individually pouched	Bulk 20 Tests/Box
8.04.18.0.0001 8.04.18.0.0100	FSH Test Strip, Urine, Individually pouched	Bulk 100 Tests/Box
8.04.19.0.0001 8.04.19.0.0015	FSH Midstream Test, Individually Pouched	Bulk 15 Tests/Box

# KIDNEY FUNCTION RAPID TESTS

## Overview

Atlas Medical Microalbumin Rapid Test is a rapid visual immunoassay used for the qualitative detection of microalbumin in human urine samples. This kit is intended for use as an aid in the diagnosis of renal dysfunction.

## Features

- 🔄 The test comes in cassette and Strip formats.



Item Code	Item Description	Sizes
8.16.52.0.0001 8.16.52.0.0020	Microalbumin Test Cassette ,Individually Pouched	Bulk 20 Tests/Box
8.16.53.0.0001 8.16.53.0.0100	Microalbumin Test Strip ,Individually Pouched	Bulk 100 Tests/Box

# INFLAMMATION AND CANCER MARKERS



## Overview

All the tests in this group are qualitative and based on lateral flow immunoassay for the detection of various inflammation and cancer markers. Some of these tests can be read by Atlas Medical rapid test reader "Quantas" to obtain quantitative results in less than 15 minutes.

Item Code	Item Description	Sizes
8.04.38.0.0020	Fecal Occult Blood Test (FOB) Test Cassette, Stool Sample, Individually Pouched	20 Tests/Box
8.04.85.0.0050	Fecal Occult Blood Test (FOB) Test Strip Stool Sample, Individually Pouched	50 Tests/Box
8.04.109.0.0020	Procalcitonin Test Cassette (PCT)	20 Tests/Box
8.16.78.0.0025	Calprotectin Test Cassette	25 Tests/Box



## Features

- Atlas Medical inflammation and cancer markers rapid tests are supplied in both cassette and strip formats.
- The kits are conveniently packed in different kit sizes of 20, 30 and 100 tests per kit.

Item Code	Item Description	Sizes
8.16.28.0.0001 8.16.28.0.0020	PSA Test Cassette, Whole Blood/Serum/Plasma), Individually Pouched	Bulk 20 Tests/Box
8.04.39.0.0001 8.04.39.0.0020	PSA Test Cassette, Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.40.0.0001 8.04.40.0.0100	PSA Test Strip, Serum/Plasma, Individually Pouched	Bulk 100 Tests/Box



# CARDIAC MARKERS RAPID TESTS

## Overview

Atlas Medical offers lateral flow immunoassay rapid tests to detect the three major cardiac markers namely: Troponin I, Myoglobin and CK-MB, as an aid in the diagnosis of myocardial infarction (MI).



## Features

- They can be used on whole blood (in addition to serum/plasma) making them ideal for emergency rooms.
- They come in single test or triple combo test cassette formats.
- The kits are conveniently packed in different kit sizes of 20, 30 tests per kit.

Item Code	Item Description	Sizes
8.04.45.0.0001 8.04.45.0.0020	Troponin I Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.46.0.0001 8.04.46.0.0020	Myoglobin Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.47.0.0001 8.04.47.0.0020	CK-MB Test Cassette, Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box
8.04.48.0.0001 8.04.48.0.0020	Cardiac Triple Test Cassette (Troponin I, CK-MB, Myoglobin), Whole Blood/Serum/Plasma, Individually Pouched	Bulk 20 Tests/Box



# DOA RAPID TESTS

## Overview

All the tests in this group are qualitative and based on lateral flow immunoassay for the detection of various inflammation and cancer markers. Some of these tests can be read by Atlas Medical rapid test reader "Quantas" to obtain quantitative results in less than 15 minutes.



## Features

- Atlas Medical inflammation and cancer markers rapid tests are supplied in both cassette and strip formats.
- The kits are conveniently packed in different kit sizes of 20, 30 and 100 tests per kit.



### STRIP AND CASSETTE FORMAT


Item Code	Item Description	Sizes
8.04.49.0.0001 8.04.49.0.0020	Morphine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.50.0.0001 8.04.50.0.0100	Morphine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.51.0.0001 8.04.51.0.0020	Marijuana (THC) Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.52.0.0001 8.04.52.0.0100	Marijuana (THC) Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.53.0.0001 8.04.53.0.0020	Amphetamine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.54.0.0001 8.04.54.0.0100	Amphetamine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.55.0.0001 8.04.55.0.0020	Barbiturates Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.56.0.0001 8.04.56.0.0100	Barbiturates Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.57.0.0001 8.04.57.0.0020	Benzodiazepines Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.58.0.0001 8.04.58.0.0100	Benzodiazepines Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.59.0.0001 8.04.59.0.0020	Cocaine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.60.0.0001 8.04.60.0.0100	Cocaine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.61.0.0001 8.04.61.0.0020	Methamphetamine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.62.0.0001 8.04.62.0.0100	Methamphetamine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.63.0.0001 8.04.63.0.0020	Methadone Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.64.0.0001 8.04.64.0.0100	Methadone Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.65.0.0001 8.04.65.0.0020	Phencyclidine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.66.0.0001 8.04.66.0.0100	Phencyclidine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box

Item Code	Item Description	Sizes
8.04.67.0.0001 8.04.67.0.0020	Tricyclic Anti-Depressants Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.04.68.0.0001 8.04.68.0.0100	Tricyclic Anti-Depressants Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.04.99.0.0001 8.04.99.0.0020	Buprenorphine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.23.0.0001 8.16.23.0.0100	Buprenorphine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.68.0.0001 8.16.68.0.0020	Tramadol Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.44.0.0001 8.16.44.0.0100	Tramadol Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.15.0.0001 8.16.15.0.0020	Methylenedioxyamphetamine (MDMA) Ecstasy Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.05.0.0001 8.16.05.0.0100	Methylenedioxyamphetamine (MDMA) Ecstasy Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.06.0.0001 8.16.06.0.0020	Opiates Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.07.0.0001 8.16.07.0.0100	Opiates Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.58.0.0001 8.16.58.0.0020	Cotinine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.59.0.0001 8.16.59.0.0100	Cotinine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.62.0.0001 8.16.62.0.0020	Oxycodone Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.63.0.0001 8.16.63.0.0100	Oxycodone Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.64.0.0001 8.16.64.0.0020	Ketamine Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.65.0.0001 8.16.65.0.0100	Ketamine Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.66.0.0001 8.16.66.0.0020	Proxyphene Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.67.0.0001 8.16.67.0.0100	Proxyphene Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.69.0.0001 8.16.69.0.0020	EDDP Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.70.0.0001 8.16.70.0.0100	EDDP Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box
8.16.60.0.0001 8.16.60.0.0020	Dolantin Test Cassette, Urine, Individually Pouched	Bulk 20 Tests/Box
8.16.61.0.0001 8.16.61.0.0100	Dolantin Test Strip, Urine, Individually Pouched	Bulk 100 Tests/Box

# DOA RAPID TESTS



Item Code	Item Description	Sizes
8.04.93.0.0001 8.04.93.0.0025	DOA Panel: 2 Drugs (Combination of any 2 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.94.0.0001 8.04.94.0.0025	DOA Panel: 3 Drugs (Combination of any 3 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.95.0.0001 8.04.95.0.0025	DOA Panel: 4 Drugs (Combination of any 4 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.79.0.0001 8.04.79.0.0025	DOA Panel: 5 Drugs (Combination of any 5 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.80.0.0001 8.04.80.0.0025	DOA Panel: 6 Drugs (Combination of any 6 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.81.0.0001 8.04.81.0.0025	DOA Panel: 7 Drugs (Combination of any 7 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.82.0.0001 8.04.82.0.0025	DOA Panel: 8 Drugs (Combination of any 8 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.83.0.0001 8.04.83.0.0025	DOA Panel: 9 Drugs (Combination of any 9 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.04.84.0.0001 8.04.84.0.0025	DOA Panel: 10 Drugs (Combination of any 10 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box
8.16.03.0.0001 8.16.03.0.0025	DOA Panel: 11 Drugs (Combination of any 11 drugs),Urine, Individually Pouched	Bulk 25 Tests/Box
8.16.04.0.0001 8.16.04.0.0025	DOA Panel: 12 Drugs (Combination of any 12 drugs), Urine, Individually Pouched	Bulk 25 Tests/Box



Item Code	Item Description	Sizes
8.42.43.0.0001	Drug Of Abuse Cup, 7 parameters (Combination of any 7 drugs), Urine, Individually Pouched	Bulk
8.42.50.0.0001	Drug Of Abuse Cup, 8 parameters (Combination of any 7 drugs), Urine, Individually Pouched	Bulk
8.16.71.0.0001	Drug Of Abuse Cup, 10 parameters (Combination of any 7 drugs), Urine, Individually Pouched	Bulk



## Quantas<sup>+</sup>™

- Quantas<sup>+</sup>™ reader is made with German technology for reliable, fast and accurate results.
- A printer can be attached to document the results.
- Can be adjusted to read both Atlas Medical Rapid tests and DOA tests.

# CLINICAL CHEMISTRY KITS

## Overview

All the tests in this group are qualitative and based on lateral flow immunoassay for the detection of various inflammation and cancer markers. Some of these tests can be read by Atlas Medical rapid test reader "Quantas" to obtain quantitative results in less than 15 minutes.

## Features

- Atlas Medical inflammation and cancer markers rapid tests are supplied in both cassette and strip formats.
- The kits are conveniently packed in different kit sizes of 20, 30 and 100 tests per kit.



Item Code	Item Description	Sizes
8.05.00.0.0250 8.05.00.0.0500	Albumin Bromocresol Green	2x125ml 4x125ml
8.05.01.0.0030 8.05.01.0.0060	Amylase	3x10ml 6x10ml
8.05.02.0.0020 8.05.02.0.0040	Acid Phosphatase Kinetic, Hilmann Method (Tablets)	10x2ml 20x2ml
8.05.03.0.0030 8.05.03.0.0060	Alkaline Phosphatase Kinetic, DGKC Method (Tablets)	10x3ml 20x3ml
8.05.04.0.0250 8.05.04.0.0500	Alkaline Phosphatase Kinetic, DGKC Method (Liquid)	5x50ml 5x100ml
8.05.05.0.0250 8.05.05.0.0500	Bilirubin Total (DMSO Method)	2x125ml 4x125ml
8.05.06.0.0250 8.05.06.0.0500	Bilirubin Direct (DMSO Method)	2x125ml 4x125ml
8.05.07.0.0250 8.05.07.0.0500	Bilirubin Total & Direct (DMSO Method)	2x125ml 4x125ml
8.05.08.0.0250 8.05.08.0.0500	Calcium Arsenazo III	2x125ml 4x125ml
8.05.09.0.0250 8.05.09.0.0500	Calcium O-Cresolphthalein	2x125ml 4x125ml
8.05.10.0.0250 8.05.10.0.0500	Chloride Thiocyanate Colorimetric	2x125ml 4x125ml
8.05.11.0.0250 8.05.11.0.0500	Cholesterol Liquid (CHOD-POD)	2x125ml 4x125ml
8.05.12.0.0025 8.05.12.0.0050	CK-MB Kinetic (Tablets)	10x2.5ml 20x2.5ml
8.05.13.0.0050 8.05.13.0.0100	CK-MB Kinetic (Liquid)	5x10ml 5x20ml
8.05.14.0.0025 8.05.14.0.0050	CK-NAC Kinetic (Tablets)	10x2.5ml 20x2.5ml
8.05.15.0.0050 8.05.15.0.0100	CK-NAC Kinetic (Liquid)	5x10ml 5x20ml



Item Code	Item Description	Sizes
8.05.16.0.0250 8.05.16.0.0500	Creatinine Jaffe Color-Kinetic	2x125ml 4x125ml
8.05.17.0.0250 8.05.17.0.0500	Glucose GOD-POD (Liquid)	2x125ml 4x125ml
8.05.18.0.0020 8.05.18.0.0040	GOT (AST) IFCC Kinetic (Tablets)	10x2ml 20x2ml
8.05.19.0.0250 8.05.19.0.0500	GOT (AST) IFCC Kinetic (Liquid)	5x50ml 5x100ml
8.05.21.0.0020 8.05.21.0.0040	GPT (ALT) IFCC Kinetic (Tablets)	10x2ml 20x2ml
8.05.22.0.0250 8.05.22.0.0500	GPT (ALT) IFCC Kinetic (Liquid)	5x50ml 5x100ml
8.05.24.0.0020 8.05.24.0.0040	Gamma GT Kinetic, Carboxy Substrate (Tablets)	10x2ml 20x2ml
8.05.25.0.0250 8.05.25.0.0500	Gamma GT Kinetic, Carboxy Substrate (Liquid)	5x50ml 5x100ml
8.05.26.0.0100 8.05.26.0.0200	HDL Cholesterol Precipitating Reagent	2x50ml 2x100ml
8.05.27.0.0250 8.05.27.0.0500	Iron Ferrozine Colorimetric	2x125ml 4x125ml
8.05.28.0.0030 8.05.28.0.0060	LDH IFCC Kinetic (Tablets)	10x3ml 20x3ml
8.05.29.0.0250 8.05.29.0.0500	LDH IFCC Kinetic (Liquid)	5x50ml 5x100ml
8.05.30.0.0060	Lipase Kinetic (Liquid)	6x10ml
8.05.31.0.0250 8.05.31.0.0500	Magnesium Calmagite Colorimetric	2x125ml 4x125ml
8.05.32.0.0250 8.05.32.0.0500	Phosphorus Phosphomolybdate UV	2x125ml 4x125ml
8.05.33.0.0050 8.05.33.0.0100	Potassium Colorimetric	50 Tests 100 Tests
8.05.34.0.0050 8.05.34.0.0100	Sodium Colorimetric	50 Tests 100 Tests
8.05.35.0.0100	TIBC (Total Iron Binding Capacity)	100 Tests
8.05.36.0.0250 8.05.36.0.0500	Total Lipids Phosphovanilline Colorimetric	2x125ml 4x125ml
8.05.37.0.0250 8.05.37.0.0500	Total Protein Biuret Colorimetric	2x125ml 4x125ml
8.05.38.0.0250 8.05.38.0.0500	Total Protein in CSF	2x125ml 4x125ml

# CLINICAL CHEMISTRY KITS



Item Code	Item Description	Sizes
8.05.39.0.0250 8.05.39.0.0500	Triglycerides GPO-POD Colorimetric	2x125ml 4x125ml
8.05.40.0.0250 8.05.40.0.0500	Urea Urease-GLDH Kinetic (UV)	5x50ml 5x100ml
8.05.41.0.0250 8.05.41.0.0500	Urea Berthelot Colorimetric	2x125ml 4x125ml
8.05.42.0.0250 8.05.42.0.0500	Uric Acid Uricase-PAP Colorimetric (Two Reagents)	2x125ml 4x125ml
8.05.45.0.0250 8.05.45.0.0500	G6PD Deficiency Qualitative Kit	250 Tests 500 Tests
8.05.45.1.0250 8.05.45.1.0500	G6PD Deficiency Qualitative Kit, (with Filter Cards)	250 Tests 500 Tests
8.05.46.0.0075 8.05.46.0.0150	G6PD Deficiency Quantitative Kit	75 Tests 150 Tests
8.05.46.1.0075 8.05.46.1.0150	G6PD Deficiency Quantitative Kit, (with Filter Cards)	75 Tests 150 Tests



Item Code	Item Description	Sizes
8.05.71.0.0250 8.05.71.0.0500	Uric Acid Uricase-PAP Colorimetric (Mono Reagents)	2x125ml 4x125ml
8.05.43.1.0005	Pathological Control for Clinical Chemistry, Lyophilized, Human Source	5ml/vial
8.05.44.1.0005	Normal Control for Clinical Chemistry, Lyophilized, Human Source	5ml/vial
8.05.47.0.0003	G6PD Control, Normal Level, (Lyophilized)	6x0.5ml
8.05.51.0.0100	HDL Cholesterol, Enzymatic Colorimetric Direct Method	100 Tests
8.05.52.0.0100	LDL Cholesterol, Enzymatic Colorimetric Direct Method	100 Tests
8.05.72.0.0250	Alkaline Phosphatase Colorimetric Method	5x50ml
8.40.00.1.0100	HbA1c Direct Enzymatic Colorimetric Kit	100 Tests
8.05.73.0.0020	Electrolytes Kit (Sodium, Potassium, Chloride)	20 Tests each
8.05.73.0.0050	Electrolytes Kit (Sodium, Potassium, Chloride)	50 Tests each

# ALCOHOL TESTS

## Overview

Atlas Medical supplies kits to test for alcohol in urine, saliva and breath. The urine alcohol test is based on the detection of EtG (Ethyl Glucuronide) in urine using a rapid lateral flow immunoassay. Whereas the saliva alcohol test uses a strip with dry chemistry pad that changes color to indicate the level of alcohol in the saliva. The alcohol breath test is a tube with crystals that change color as the subject blows through when alcohol level in breath exceeds a certain limit.



Item Code	Item Description	Sizes
8.25.01.0.0001 8.25.01.0.0100	Saliva Alcohol Test Strip, Individually Pouched	Bulk 100 Tests/Box
8.25.02.0.0001 8.25.02.0.0025	Breath Alcohol Test (0.02%), Individually Pouched	Bulk 25 Tests/Box
8.25.03.0.0001 8.25.03.0.0025	Breath Alcohol Test (0.05%), Individually Pouched	Bulk 25 Tests/Box
8.25.04.0.0001 8.25.04.0.0025	Breath Alcohol Test (0.08%), Individually Pouched	Bulk 25 Tests/Box
8.25.05.0.0001 8.25.05.0.0025	Ethyl glucuronide (EtG) Test Strip, Individually Pouched	Bulk 25 Tests/Box

## Features

- The kits are conveniently packed in different kit sizes of 25, 50 and 100 tests per kit.



# STAINS FOR HISTOLOGY & MICROBIOLOGY



## Overview

Atlas Medical is well known for its range of lab stains for histology and microbiology applications.

Atlas Medical stains are made of the highest quality ingredients to ensure good quality and vivid staining.

## Features

- The stains come in convenient sizes, but custom sizes are also available.



### STAIN PACKS FOR HISTOLOGY

Item Code	Item Description	Sizes
8.17.003.0300	Periodic Acid Schiff (PAS) Stain Kit	3x100ml
8.17.004.0300	Iron Stain Kit - Perl	3x100ml
8.17.009.1000	Gram Stain Pack	4x250ml
8.17.010.0750	Cold ZN - Kinyoun Stain Pack	3x250ml
8.17.011.0750	ZN Pack Standard	3x250ml
8.17.015.0500	Diff-3 Stain Pack	4x125ml
8.17.016.1000	Papanicolaou Stain Pack (EA35, EA50, EA65, OG6)	4X250ml

### STAINS FOR HISTOLOGY

Item Code	Item Description	Sizes
8.15.017.0250	Carbol Fuchsin (Gram)	250ml/Bottle
8.15.019.0250	Carbol Fuchsin (Ziehl-Neelsen)	250ml/Bottle
8.15.032.0250	Crystal Violet (for Gram Stain)	250ml/Bottle
8.15.037.0250	Eosin Y (1% Aqueous)	250ml/Bottle
8.15.038.0250	Eosin Y (5% Aqueous)	250ml/Bottle
8.15.039.0250	Eosin Y (1% Alcoholic)	250ml/Bottle
8.15.041.0250	Field Stain (Solution A)	250ml/Bottle
8.15.042.0250	Field Stain (Solution B)	250ml/Bottle
8.15.047.0250	Giemsa Stain (Modified-Glycerol/Methanol)	250ml/Bottle
8.15.049.0250	Gram's Iodine	250ml/Bottle
8.15.051.0250	Gram's Decolouriser	250ml/Bottle
8.15.059.0250	Haematoxylin Harris (no Acetic Acid)	250ml/Bottle
8.15.060.0250	Haematoxylin Harris (with Acetic Acid)	250ml/Bottle
8.15.069.0250	Leishman Stain	250ml/Bottle
8.15.074.0250	Lugol's Iodine	250ml/Bottle
8.15.076.0250	Malachite Green (Aqueous)	250ml/Bottle
8.15.078.0250	May Grunwald Stain (Modified)	250ml/Bottle
8.15.105.0250	New Methylene Blue for Reticulocytes	250ml/Bottle
8.15.110.0250	Papanicolaou Stain EA35	250ml/Bottle
8.15.111.0250	Papanicolaou Stain EA36	250ml/Bottle
8.15.112.0250	Papanicolaou Stain EA65	250ml/Bottle
8.15.114.0250	Papanicolaou Stain EA50	250ml/Bottle
8.15.115.0250	Papanicolaou Stain OG6	250ml/Bottle
8.15.126.0250	Safranin (1% Aqueous)	250ml/Bottle
8.15.143.0250	Wright's Stain (Modified)	250ml/Bottle
8.15.144.0250	ZN Decolouriser	250ml/Bottle
8.15.146.0100	Immersion Oil	100ml/Bottle



### MICROBIOLOGY

Item Code	Item Description	Sizes
8.38.00.0.0025	Blood Culture Bottles, Pediatric Size	25ml/Bottle
8.38.00.0.0050	Blood Culture Bottles, Adult Size	50ml/Bottle
8.36.00.0.0020	Mycoplasma IES Infection Test	20 Tests/Box
8.36.01.0.0030	Mycoplasma IES Plus Infection Test	30 Tests/Box

# ANTIBIOTIC SENSITIVITY DISCS

## Overview

Antibiotic sensitivity is a term used to describe the susceptibility of bacteria to antibiotics. Antibiotic susceptibility testing (AST) is usually carried out to determine which antibiotic will be most successful in treating a bacterial infection in vivo.

Small discs containing antibiotics are placed onto a plate upon which bacteria are growing. If the bacteria are sensitive to the antibiotic, a clear ring, or zone of inhibition, is seen around the disc indicating poor growth.

## Features

- Atlas Medical offers a wide range of antibiotics discs at competitive prices.
- Easy to use.
- The kit comes with a Cartridge Applicator.
- Reliable quality.
- Comprehensive range of antibiotics at different concentrations.



Item Code	Item Description	Sizes
8.39.01.0.0250	AMIKACIN (30 µg) – AK	5x50 Discs
8.39.02.0.0250	AMOXICILLIN (10 µg) – AX	5x50 Discs
8.39.03.0.0250	AMOXICILLIN / CLAVULANIC ACID (20 µg + 10 µg) - AC	5x50 Discs
8.39.04.0.0250	AMPICILLIN (10 µg) – AP	5x50 Discs
8.39.05.0.0250	AMPICILLIN / SULBACTAM (10 µg - 10 µg) – AS	5x50 Discs
8.39.06.0.0250	AZITHROMYCIN (15 µg) – AZ	5x50 Discs
8.39.07.0.0250	AZTREONAM (30 µg) – AT	5x50 Discs
8.39.08.0.0250	CEFACLOR (30 µg) - CG	5x50 Discs
8.39.09.0.0250	CEFADROXIL (30 µg) - CD	5x50 Discs
8.39.10.0.0250	CEFAZOLIN (30 µg) - CF	5x50 Discs
8.39.11.0.0250	CEFDINIR (5µg) - CN	5x50 Discs
8.39.12.0.0250	CEFIXIME (5 µg) - FX	5x50 Discs
8.39.13.0.0250	CEFOPERAZONE (75 µg) - PZ	5x50 Discs
8.39.14.0.0250	CEFOPERAZONE / SULBACTAM (75 µg + 30 µg) - CS	5x50 Discs
8.39.15.0.0250	CEFOTAXIME (30 µg) - CX	5x50 Discs
8.39.16.0.0250	CEFPIROME (30 µg) - CE	5x50 Discs
8.39.17.0.0250	CEFPODOXIME (10 µg) – CO	5x50 Discs
8.39.18.0.0250	CEFPROZIL (30 µg) - FP	5x50 Discs
8.39.19.0.0250	CEFTAZIDIME (30 µg) – CZ	5x50 Discs
8.39.20.0.0250	CEFTIZOXIME (30 µg) - FO	5x50 Discs
8.39.21.0.0250	CEFTRIOXONE (30 µg) - FR	5x50 Discs
8.39.22.0.0250	CEFUROXIME (30 µg) - CR	5x50 Discs
8.39.23.0.0250	CEPHALEXIN (30 µg) - CP	5x50 Discs
8.39.25.0.0250	CEPHALOTHIN (30 µg) - CA	5x50 Discs
8.39.26.0.0250	CHLORAMPHENICOL (30 µg) - CK	5x50 Discs
8.39.27.0.0250	CIPROFLOXACIN (5 µg) - CL	5x50 Discs
8.39.28.0.0250	CLARITHROMYCIN (15 µg) - CL	5x50 Discs
8.39.29.0.0250	CLINDAMYCIN (2 µg) - CM	5x50 Discs
8.39.30.0.0250	CLOXACILLIN (5 µg) - CV	5x50 Discs
8.39.32.0.0250	DOXYCYCLINE (30 µg) - DO	5x50 Discs
8.39.33.0.0250	ERYTHROMYCIN (15 µg) - ER	5x50 Discs
8.39.34.0.0250	FURAZOLIDONE (100 µg) - FZ	5x50 Discs
8.39.35.0.0250	GATIFLOXACIN (5 µg) - GF	5x50 Discs
8.39.36.0.0250	GENTAMYCIN (10 µg) - GM	5x50 Discs
8.39.38.0.0250	KANAMYCIN (30 µg) - KA	5x50 Discs
8.39.39.0.0250	LEVOFOLXACIN (5 µg) - LV	5x50 Discs
8.39.40.0.0250	LINCOMYCIN (15 µg) - LN	5x50 Discs
8.39.41.0.0250	LINEZOLID (30 µg) - LI	5x50 Discs
8.39.42.0.0250	LOMEFLOXACIN (10 µg) - LF	5x50 Discs
8.39.43.0.0250	MEROPENEM (10 µg) - MR	5x50 Discs
8.39.44.0.0250	MINOCYCLINE (30 µg) - MN	5x50 Discs
8.39.45.0.0250	MOXIFLOXACIN (5 µg) - MF	5x50 Discs
8.39.46.0.0250	NALIDIXIC ACID (30 µg) - NA	5x50 Discs
8.39.47.0.0250	NITROFURANTOIN (300 µg) - FU	5x50 Discs

Item Code	Item Description	Sizes
8.39.48.0.0250	NORFLOXACIN (10 µg) - NF	5x50 Discs
8.39.49.0.0250	OFLOXACIN (5 µg) - OF	5x50 Discs
8.39.50.0.0250	PEFLOXACIN (5 µg) – PF	5x50 Discs
8.39.51.0.0250	PENICILLIN –G (10 IU) – PG	5x50 Discs
8.39.52.0.0250	PIPERACILLIN (100 µg) – PC	5x50 Discs
8.39.53.0.0250	PIPERACILLIN / TAZOBACTAM (100 µg + 10 µg) - PT	5x50 Discs
8.39.54.0.0250	RIFAMPIN (5 µg) – RN	5x50 Discs
8.39.55.0.0250	ROXITHROMYCIN (30 µg) – RO	5x50 Discs

# ANTIBIOTIC SENSITIVITY DISCS



Item Code	Item Description	Sizes
8.39.56.0.0250	SPARFLOXACIN (5 µg) – SP	5x50 Discs
8.39.57.0.0250	STREPTOMYCIN (10 µg) – ST	5x50 Discs
8.39.59.0.0250	TEICOPLANIN (30 µg) – TC	5x50 Discs
8.39.60.0.0250	TETRACYCLINE (30 µg) – TE	5x50 Discs
8.39.61.0.0250	TICARCILLIN / CLAVULANIC ACID (75 µg + 2.5 µg)-TC	5x50 Discs
8.39.62.0.0250	TOBRAMYCIN (10 µg) – TO	5x50 Discs
8.39.63.0.0250	TRIMETHOPRIM (5 µg) – TR	5x50 Discs
8.39.64.0.0250	VANCOMYCIN (30 µg) – VM	5x50 Discs
8.39.65.0.0250	POLYMYXIN-B (300 UNITS) -PB	5x50 Discs



# ELISA KITS

## Features

- ⌚ The kits feature high sensitivities, simple and robust methods, breakable well strips, quantitative results, ready-to use liquid reagents, and reasonable assay time.
- ⌚ The assays can be used on most open ELISA manual readers and washers as well as open ELISA auto-analyzers.
- ⌚ Kits are packed in sizes of 96 tests.



## Overview

Atlas Medical offers a range of Enzyme Linked Immunosorbent Assay (ELISA or EIA) to detect major hormones in the fields of thyroids and fertility in serum.



HORMONE ELISA KITS		
Item Code	Item Description	Sizes
8.07.02.0.0096	PSA Elisa Kit	96 Tests
8.07.10.0.0096	Free PSA Elisa Kit	96 Tests



HORMONE ELISA KITS		
Item Code	Item Description	Sizes
8.10.01.0.0096	hCG Elisa Kit	96 Tests
8.10.03.0.0096	FSH Elisa Kit	96 Tests
8.10.04.0.0096	LH Elisa Kit	96 Tests
8.10.05.0.0096	Prolactin Elisa Kit	96 Tests
8.12.00.0.0096	T3 Elisa Kit	96 Tests
8.12.01.0.0096	T4 Elisa Kit	96 Tests
8.12.02.0.0096	TSH Elisa Kit	96 Tests
8.12.03.0.0096	Free T4 Elisa Kit	96 Tests
8.12.04.0.0096	Free T3 Elisa Kit	96 Tests
8.07.09.0.0096	Free β-hCG Elisa Kit	96 Tests
8.07.09.1.0096	Beta hCG Elisa Kit	96 Tests
8.11.03.0.0096	Progesterone Elisa kit	96 Tests
8.11.04.0.0096	Testosterone Elisa Kit	96 Tests

# ELISA KITS

## Overview

Atlas Medical offers a range of Enzyme Linked Immunosorbent Assay (ELISA or EIA) to detect IgG and IgM antibodies against ToRCH (Toxoplasmosis, Rubella, CMV and Herpes I & II) in serum.

- \* Different Packaging sizes.
- \* Easy to Use
- \* High Quality



Item Code	Item Description	Sizes
8.13.10.0.0096	Herpes Simplex 2 IgM (HSV2 IgM) Elisa Kit	96 Tests
8.13.11.0.0096	Herpes Simplex 1,2 IgG (HSV1,2 IgG) Elisa Kit	96 Tests
8.13.12.0.0096	Herpes Simplex 1,2 IgM (HSV1,2 IgM) Elisa Kit	96 Tests



### TORCH ELISA KITS

Item Code	Item Description	Sizes
8.13.00.0.0096	Toxo Plasma Gondii IgG (Toxo IgG) Elisa kit	96 Tests
8.13.01.0.0096	Toxoplasma gondii IgM (Toxo IgM) Elisa kit	96 Tests
8.13.02.0.0096	Rubella IgG Elisa Kit	96 Tests
8.13.03.0.0096	Rubella IgM Elisa Kit	96 Tests
8.13.05.0.0096	Cytomegalovirus IgG (CMV IgG) Elisa Kit	96 Tests
8.13.06.0.0096	Cytomegalovirus IgM (CMV IgM) Elisa Kit	96 Tests
8.13.07.0.0096	Herpes Simplex 1 IgG (HSV1 IgG) Elisa Kit	96 Tests
8.13.08.0.0096	Herpes Simplex 1 IgM (HSV1 IgM) Elisa Kit	96 Tests
8.13.09.0.0096	Herpes Simplex 2 IgG (HSV2 IgG) Elisa Kit	96 Tests

## Overview

Atlas Medical offers a range of Enzyme Linked Immunosorbent Assay (ELISA or EIA) to detect a series of infection diseases such as HIV, Hepatitis (A, B, C, D and E) and H. pylori (antigens in feces).



### INFECTIOUS DISEASES ELISA KITS

Item Code	Item Description	Sizes
8.14.19.1.0096	Helicobacter Pylori Antigen Elisa Kit	96 Tests



### INFECTIOUS DISEASES ELISA KITS

Item Code	Item Description	Sizes
8.14.28.0.0096	HBsAg Elisa Kit	96 Tests
8.14.29.0.0096	HBsAb Elisa Kit	96 Tests
8.14.30.0.0096	HBCAb Elisa Kit	96 Tests
8.14.31.0.0096	HBeAg Elisa Kit	96 Tests
8.14.32.0.0096	HBeAb Elisa Kit	96 Tests
8.14.38.0.0096	HCV Ab Elisa Kit	96 Tests
8.14.39.0.0096	HAV IgM Elisa Kit	96 Tests
8.14.40.0.0096	HIV 1,2 Antibody Elisa Kit	96 Tests



### OTHER ELISA KITS

Item Code	Item Description	Sizes
8.07.03.0.0096	Alpha Feto Protein (AFP) Elisa Kit	96 Tests
8.07.08.0.0096	Ferritin Elisa Kit	96 Tests
8.08.00.0.0096	Troponin I Elisa Kit	96 Tests
8.09.00.0.0096	IgE Elisa Kit	96 Tests
8.51.00.0.0096	25-OH Vitamin D Elisa Kit	96 Tests
8.57.00.0.0096	Vitamin B12 Elisa Kit	96 Tests
8.58.00.0.0096	Folic Acid Elisa Kit	96 Tests
8.06.32.0.0096	Anti-CRA Elisa K	96 Tests



# HOME TESTS



## Overview


Atlas Medical provides a range of home tests that have been specifically CE marked for OTC use. The range includes fertility tests (Pregnancy, Ovulation and Menopause). The home tests range also includes other medical conditions such as liver function, kidney function, diabetes and urine tract infection. These tests are based on urine reagent strips.

### Screening Kits

Item Code	Item Description	Sizes
70004001	Atlas Home Diabetes Test	2 Tests/Box
70021001	Atlas Home Urinary Tract Infection Test	2 Tests/Box
70022001	Atlas Home Kidney Function Test	2 Tests/Box
70023001	Atlas Home Liver Function Test	2 Tests/Box

## Features

- These tests come in cassette, midstream and strip formats.
- The screening kits come with 2 individually pouched strips and easy to read instructions for use.
- All kits are packed in attractively designed boxes with various languages.
- Atlas Medical also supplies these kits under OEM arrangements.



- Simply dip the test strip in urine for one second.
- Wait for 30 - 60 seconds.
- Compare the colors on the strip to the color chart on the box.



### Fertility Kits

Item Code	Item Description	Sizes
70171001	Atlas Home Pregnancy Test Cassette	1 Test/Box
70172001	Atlas Home Pregnancy Test Midstream	1 Tests/Box
70174001	Atlas Home Ovulation Test Cassette	5 Tests/Box
70175001	Atlas Home Ovulation Test Midstream	3 Tests/Box
70177001	Atlas Home Menopause Test Cassette	1 Test/Box
70178001	Atlas Home Menopause Test Midstream	1 Test/Box
70180001	Atlas Home Pregnancy Test Strip (With Handle)	1 Test/Box
70170001	Atlas Home Pregnancy Test Strip	1 Test/Box

# BLOOD GLUCOSE MONITORING SYSTEMS

## Overview

Testing your blood glucose regularly helps you better manage your diabetes. Reliance™ by Atlas Medical, uses the latest blood glucose sensor technologies to offer you the most accurate and reliable results for the peace of mind you need. Atlas Medical offers these systems in different Electrode strips types which includes Gold and Silver Electrodes.



## Features

- Reliance Gold™ by Atlas Medical, uses the latest blood glucose gold sensor technology to offer the most accurate and reliable results.
- Test time required is 5 Seconds.
- Required sample volume is 0.9µl
- Test result range is between 10 - 600 mg/dl (0.6 -33.3 mmol/L)

### Reliance Gold

Item Code	Item Description	Sizes
8.52.00.0.0001	Reliance Gold Glucometer Pack	1 Pack
8.52.00.0.0050	Strips for Reliance Gold Glucometer	50 Strips/Bottle

## Features

- Test time required is 5 Seconds.
- Required sample volume is 0.6µl
- Test result range is between 20 - 600 mg/dl (1.3 - 33.3 mmol/L)



### Reliance Silver

Item Code	Item Description	Sizes
8.53.00.0.0001	Reliance Silver Glucometer Pack	1 Pack
8.53.00.0.0050	Strips for Reliance Silver Glucometer	50 Strips/Bottle



# CERTIFICATES



ISO 13485



Blood Grouping CE Certificate



Full Quality Assurance Certificate

## OTHER CERTIFICATES

- FDA 510k Atlas Drug of Abuse Tests ( Cup & Panel Format )
- GMP Certificate
- FDA 510k Atlas Home Pregnancy Test (Midstream Format)
- hCG Test Strip CE certificate
- hCG Test Cassette CE certificate
- hCG Midstream Test CE certificate
- Ovulation Test Midstream CE certificate
- Ovulation Test Cassette CE certificate
- Menopause Test Midstream CE certificate
- Menopause Test Cassette CE certificate
- Liver Function Test CE certificate
- Diabetes Test CE certificate
- UTI Test CE certificate
- Kidney function Test CE certificate

# CONTACT US

---



## FIND US ON:

---



@Atlas.Medical.UK



[linkedin.com/company/atlas-medical-co](https://www.linkedin.com/company/atlas-medical-co)



@AtlasMedicalUK



[Amazon.co.uk/Atlas Medical](https://www.amazon.co.uk/Atlas-Medical)

# INTERNATIONAL PRESENCE

---

