

ComASP® Cefiderocol 0.008-128 µg/mL

Broth microdilution panel for Cefiderocol susceptibility testing according to ISO 20776-1



ComASP® Cefiderocol 0.008-128 µg/mL is a 2-test panel containing the dried-up antibiotic in 15 two-fold dilutions for the quantitative determination of the minimum inhibitory concentration (MIC) of Cefiderocol against Gram-negative non-fastidious organisms such *Enterobacterales* and *Pseudomonas aeruginosa*, *Acinetobacter* spp., *Stenotrophomonas maltophilia*.

Kit content:

- 4 individually packed ComASP® Cefiderocol 0.008-128 µg/mL panels (2 samples/panel)
- 8 Tubes of **Iron-Depleted Cation-adjusted Mueller Hinton Broth (ID-CAMHB)** (3.6 ml)
- Sealing Film

Ref.	Device format	Description	Packaging
75009	ComASP®	Cefiderocol 0.008-128 µg/mL	8 tests




Liofilchem srl
Headquarters, Manufacturing Site,
International distribution
Roseto degli Abruzzi, Italy



Liofilchem, Inc.
US Distribution Center
Waltham, MA, USA



Liofilchem Nordics ApS
Scandinavia Customer Center
Copenhagen, Denmark



MTS™ Cefiderocol
0.016-256 µg/mL
for P. aeruginosa



Cefiderocol disc
30 µg

Ref.	Device format	Description	Packaging
92067	MTS™ (MIC Test Strip)	Cefiderocol 0.016-256 µg/mL	30 mic test
920670	MTS™ (MIC Test Strip)	Cefiderocol 0.016-256 µg/mL	100 mic test
920671	MTS™ (MIC Test Strip)	Cefiderocol 0.016-256 µg/mL	10 mic test
9266	Antibiotic disc in cartridge	Cefiderocol 30 µg	5x50 Discs
9266/1	Antibiotic disc in cartridge	Cefiderocol 30 µg	50 Discs
9266/2	Antibiotic disc in canister	Cefiderocol 30 µg	250 Discs



IVD



Liofilchem srl
 Headquarters, Manufacturing Site,
 International distribution
 Roseto degli Abruzzi, Italy



Liofilchem, Inc.
 US Distribution Center
 Waltham, MA, USA



Liofilchem Nordics ApS
 Scandinavia Customer Center
 Copenhagen, Denmark