

## ANTIBIOTIC DISCS IN CARTRIDGES

DESCRIPTION		µg	CLSI 1,2	EUCAST 3	BSAC 4	PACKAGING	REF.
-------------	--	----	----------	----------	--------	-----------	------



Amikacin	AK	30	✓	✓	✓	5x50 Discs	9004
Amoxicillin	AML	30				5x50 Discs	9005
Amoxicillin	AML	10			✓	5x50 Discs	9133
Amoxicillin	AML	2		✓	✓	5x50 Discs	9151
Ampicillin	AMP	2		✓	✓	5x50 Discs	9115
Ampicillin	AMP	10	✓	✓	✓	5x50 Discs	9006
Ampicillin / Sulbactam	AMS	20 (10+10)	✓	✓		5x50 Discs	9031
Ampiclox (AMP / CX)	ACL	30 (25+5)				5x50 Discs	9122
Amoxicillin/Clavulanic acid	AUG	30 (20+10)	✓	✓	✓	5x50 Discs	9048
Azithromycin	AZM	15	✓	✓	✓	5x50 Discs	9105
Azlocillin	AZL	75	✓			5x50 Discs	9007
Aztreonam	ATM	30	✓	✓	✓	5x50 Discs	9008
Bacitracin	BA	10 IU				5x50 Discs	9051
Carbenicillin	CAR	100	✓		✓	5x50 Discs	9009
Cefaclor	CEC	30	✓	✓	✓	5x50 Discs	9010
Cefadroxil	CDX	30		✓	✓	5x50 Discs	9052
Cefamandole	MA	30	✓		✓	5x50 Discs	9014
Cefazolin	KZ	30	✓	✓		5x50 Discs	9015
Cefepime	FEP	30	✓	✓	✓	5x50 Discs	9104
Cefixime	CFM	5	✓	✓	✓	5x50 Discs	9089
Cefoperazone	CFP	75	✓			5x50 Discs	9108
Cefoperazone	CFP	30			✓	5x50 Discs	9016
Cefotaxime	CTX	30	✓		✓	5x50 Discs	9017
Cefotaxime	CTX	5		✓	✓	5x50 Discs	9152
Cefotetan	CTT	30	✓			5x50 Discs	9081
Cefpodoxime	PX	10	✓	✓	✓	5x50 Discs	9064
Cefoxitin	FOX	30	✓	✓	✓	5x50 Discs	9018
Cefprozil	CPR	30	✓			5x50 Discs	9112
Cefsulodin	CSD	30				5x50 Discs	9053
Ceftazidime	CAZ	30	✓		✓	5x50 Discs	9019
Ceftazidime	CAZ	10		✓		5x50 Discs	9153
Ceftibuten	CTB	30	✓	✓		5x50 Discs	9101

## ANTIBIOTIC DISCS IN CARTRIDGES

DESCRIPTION		µg	CLSI <sup>1,2</sup>	EUCAST <sup>3</sup>	BSAC <sup>4</sup>	PACKAGING	REF.
Ceftizoxime	CZX	30	✓		✓	5x50 Discs	9054
Ceftriaxone	CRO	30	✓	✓	✓	5x50 Discs	9020
Cefuroxime	CXM	30	✓	✓	✓	5x50 Discs	9021
Cephalexin	CL	30			✓	5x50 Discs	9011
Cephalothin	KF	30	✓		✓	5x50 Discs	9013
Cefazolin	KZ	30	✓			5x50 Discs	9015
Cephradine	CE	30			✓	5x50 Discs	9055
Chloramphenicol	C	30	✓	✓	✓	5x50 Discs	9022
Chloramphenicol	C	10			✓	5x50 Discs	9128
Cinoxacin	CIN	100	✓			5x50 Discs	9057
Ciprofloxacin	CIP	5	✓	✓	✓	5x50 Discs	9056
Clarithromycin	CLR	15	✓	✓		5x50 Discs	9098
Clindamycin	CD	2	✓	✓	✓	5x50 Discs	9047
Cloxacillin	CX	5				5x50 Discs	9058
Colistin sulfate	CS	10	✓			5x50 Discs	9023
Dicloxacillin	DCX	1				5x50 Discs	9093
Doripenem	DOR	10	✓	✓		5x50 Discs	9154
Doxycycline	DXT	30	✓		✓	5x50 Discs	9059
Ertapenem	ETP	10	✓	✓	✓	5x50 Discs	9061
Erythromycin	E	15	✓	✓		5x50 Discs	9024
Fosfomicin (includes G-6-P)	FOS	200	✓			5x50 Discs	9109
Fosfomicin (includes G-6-P)	FOS	100				5x50 Discs	9121
Fosfomicin	FOS	50			✓	5x50 Discs	9025
Furazolidon	FR	50				5x50 Discs	9099
Fusidic acid	FC	10		✓	✓	5x50 Discs	9049
Fusidic acid	FC	30				5x50 Discs	9111
Gatifloxacin	GAT	5	✓			5x50 Discs	9169
Gentamicin	CN	10	✓	✓	✓	5x50 Discs	9026
Gentamicin	CN	30		✓		5x50 Discs	9125
Gentamicin	CN	120				5x50 Discs	9124
Imipenem	IMI	10	✓	✓	✓	5x50 Discs	9079
Imipenem + Cilastatin	IMC	20 (10+10)				5x50 Discs	9095
Kanamycin	K	30	✓			5x50 Discs	9027
Levofloxacin	LEV	5	✓	✓	✓	5x50 Discs	9102
Lincomycin	MY	15				5x50 Discs	9116
Lincomycin	MY	2				5x50 Discs	9028
Linezolid	LNZ	10		✓	✓	5x50 Discs	9155
Linezolid	LNZ	30	✓			5x50 Discs	9136
Lomefloxacin	LOM	10	✓			5x50 Discs	9113
Loracarbef	LOR	30	✓			5x50 Discs	9173
Mecillinam	MEC	10	✓	✓	✓	5x50 Discs	9156
Meropenem	MRP	10	✓	✓	✓	5x50 Discs	9068
Metronidazole	LZ	50				5x50 Discs	9119
Metronidazole	LZ	5			✓	5x50 Discs	9076
Mezlocillin	MEZ	75	✓		✓	5x50 Discs	9062
Minocycline	MN	30	✓	✓	✓	5x50 Discs	9030
Miokamicin	MK	15				5x50 Discs	9106
Moxifloxacin	MOX	5	✓	✓	✓	5x50 Discs	9103

## ANTIBIOTIC DISCS IN CARTRIDGES

DESCRIPTION		µg	CLSI <sup>1,2</sup>	EUCAST <sup>3</sup>	BSAC <sup>4</sup>	PACKAGING	REF.
Mupirocin	MUP	200		✓		5x50 Discs	9157
Nafcillin	NAF	1	✓			5x50 Discs	9174
Nalidixic acid	NA	30	✓	✓	✓	5x50 Discs	9001
Neomycin	N	30				5x50 Discs	9032
Netilmicin	NET	10		✓		5x50 Discs	9170
Netilmicin	NET	30	✓			5x50 Discs	9033
Nitrofurantoin	F	100		✓		5x50 Discs	9158
Nitrofurantoin	F	300	✓			5x50 Discs	9034
Norflloxacin	NOR	10	✓	✓		5x50 Discs	9035
Novobiocin	NO	30				5x50 Discs	9063
Novobiocin	NO	5				5x50 Discs	9117
Ofloxacin	OFX	5	✓	✓	✓	5x50 Discs	9080
Oxacillin	OX	1	✓	✓	✓	5x50 Discs	9036
Oxacillin	OX	5				5x50 Discs	9135
Oxolinic acid	OA	2				5x50 Discs	9002
Oxytetracycline	OT	30				5x50 Discs	9065
Pefloxacin	PEF	5				5x50 Discs	9091
Penicillin G	P	10 IU	✓			5x50 Discs	9037
Penicillin G	P	2 IU				5x50 Discs	9127
Penicillin G	P	1 IU		✓	✓	5x50 Discs	9130
Phenoxymethylpenicillin	PV	10		✓		5x50 Discs	9171
Pipemidic acid	PI	20				5x50 Discs	9003
Piperacillin	PRL	30		✓		5x50 Discs	9159
Piperacillin	PRL	100	✓			5x50 Discs	9038
Piperacillin/tazobactam	TZP	110 (100+10)	✓			5x50 Discs	9100
Piperacillin/tazobactam	TZP	36 (30+6)		✓		5x50 Discs	9160
Polimyxin B	PB	100 IU				5x50 Discs	9066
Polimyxin B	PB	300 IU				5x50 Discs	9120
Quinupristin/Dalfopristin	QDA	15	✓	✓	✓	5x50 Discs	9161
Rifampicin	RD	30				5x50 Discs	9039
Rifampicin	RD	5		✓	✓	5x50 Discs	9118
Roxitromycin	RXT	15				5x50 Discs	9060
Sisomicin	SIS	30				5x50 Discs	9046
Sodium Fusidate	FC	30				5x50 Discs	9131
Spectinomycin	SPC	100	✓			5x50 Discs	9067
Spiramycin	SP	100				5x50 Discs	9088
Streptomycin	S	10	✓		✓	5x50 Discs	9040
Streptomycin	S	300		✓		5x50 Discs	9162
Sulbactam	SU	20				5x50 Discs	9129
Sulfafurazole	SF	300				5x50 Discs	9041
Sulfamethoxazole	SMX	50				5x50 Discs	9084
Sulfaprim	SXT	50				5x50 Discs	9132
Sulphonamide	S3	300				5x50 Discs	9126
Teicoplanin	TEC	30	✓	✓	✓	5x50 Discs	9050
Telithromycin	TEL	15	✓	✓	✓	5x50 Discs	9172
Tetracycline	TE	30	✓	✓		5x50 Discs	9043
Tiamulin	T	30				5x50 Discs	9094
Ticarcillin	TC	75	✓	✓	✓	5x50 Discs	9070

## ANTIBIOTIC DISCS IN CARTRIDGES

DESCRIPTION		µg	CLSI 1,2	EUCAST 3	BSAC 4	PACKAGING	REF.
Ticarcillin /Clavulanic acid	TTC	85 (75+10)	✓	✓		5x50 Discs	9096
Tigecycline	TGC	15	✓	✓	✓	5x50 Discs	9147
Tobramycin	TOB	10	✓	✓	✓	5x50 Discs	9044
Tobramycin	TOB	30		✓		5x50 Discs	9163
Trimethoprim	TM	5	✓	✓	✓	5x50 Discs	9110
Trimethoprim	TM	2.5			✓	5x50 Discs	9083
Trimethoprim-Sulfamethoxazole	SXT	25 (1,25+23,75)	✓	✓	✓	5x50 Discs	9042
Tylosin	TY	30				5x50 Discs	9082
Vancomycin	VA	5		✓	✓	5x50 Discs	9164
Vancomycin	VA	30	✓			5x50 Discs	9045

Single cartridges of 50 discs are available: add /1 to the catalogue ref. no. to indicate the relevant item.

Example: ref. 9045/1 indicates Vancomycin 30 µg in one single cartridge of 50 discs.

**NOTE:**

1. CLSI M100-S20 - Performance Standards for Antimicrobial Susceptibility Testing.

2. CLSI M02-A10 - Performance Standards for Antimicrobial Disk Susceptibility Tests; Approved Standard-Tenth Edition

3. EUCAST Clinical Breakpoint Table v. 1.3 2011-01-05

4. BSAC Methods for Antimicrobial Susceptibility Testing Version 10.1 April 2011