



*Assure Tech,
Infectious Disease Diagnostic Specialist,
Your Reliable Partner*

Dengue Rapid Test



As an aid in the diagnosis of primary and secondary Dengue infections

- ✓ Stable & Reliable Products
- ⚡ Easy to use
- ⚡ Superior Performance



Assure Tech. (Hangzhou) Co., Ltd.
Building 4, No. 1418-50, Moganshan Road,
Gongshu District, Hangzhou, 310011 Zhejiang, P.R. China



400-853-5577



Dengue Rapid Test

Dengue Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of IgG and IgM antibodies to Dengue virus in human whole blood, serum, or plasma.



Notable Features


- Reliable & Inexpensive
- Easy-to-follow procedure
- Better performance makes the result more accurate
- Can be performed by untrained or minimally skilled personnel



Product Information

Product	REF	Specimen	Test Time	Test Format	Kit Size	Qualification
Dengue Rapid Test	DEN-W23	W/S/P	15 minutes	Device	25,40 pcs/box	CE

Test Procedure



Assure Tech. (Hangzhou) Co., Ltd.
 Building 4, No. 1418-50, Moganshan Road,
 Gongshu District, Hangzhou, 310011 Zhejiang, P.R. China



400-853-5577