

# H. pylori Antigen Rapid Test (Feces)

For laboratory professional use

As an aid in the diagnosis of H. pylori infection in symptomatic or asymptomatic individuals

H. pylori Antigen Rapid Test is a rapid immunochromatographic assay for the qualitative, presumptive detection of Helicobacter pylori antigens in human fecal specimens.



## Notable Features

- Easy to Use
- Fast Detection: Result in 10 minutes.
- Rapid & accurate test results
- High Accuracy.

## Product Information

Information Type	Specimen	Detection	Format	Package	Time to results	Shelf life & Storage Conditions
Product Detail	Fecal	<i>Helicobacter pylori</i> antigens	Device	25 pcs/box	10 minutes	24 months at 2-30 C

## Test Procedure

