



PREGNANCY TEST



- Rapid, Visual Chromatographic Assay for Detection of Early Pregnancy
- Sensitivity: >99.9%, Specificity: >99.9%
- Multiple Independent clinical studies
- Reading time: 3-5minutes
- Storage: 2-30°C

Product Name	Catalog No.	Specimen Used	Test Format	CE mark	Size	US FDA
Pregnancy Test	HCG-U11	Urine	Strip	CE	50 pcs/box	
	HCG-C11	Urine/Serum	Strip	CE	50 pcs/box	√FDA 510k

Test Procedure

HCG-U11




HCG-C11





PREGNANCY TEST




- Rapid, Visual Chromatographic Assay for Detection of Early Pregnancy
- Sensitivity: >99.9%, Specificity: >99.9%
- Multiple Independent clinical studies
- Reading time: 3-5minutes
- Storage: 2-30°C

Product Name	Catalog No.	Specimen Used	Test Format	CE mark	Size	US FDA
Pregnancy Test	HCG-U23	Urine	Device	CE	40 pcs/box	
	HCG-C23	Urine/Serum	Device	CE	40 pcs/box	√FDA 510k

Test Procedure

HCG-U23



HCG-C23

