



Assure Tech,
Infectious Disease
Diagnostic Specialist,
Your Reliable Partner

Rotavirus and Adenovirus Combo Rapid Test

Rapid Screening Test for Rotavirus and Adenovirus infection



Instant Screening Test



Easy & Rapid



High Performance

Assure Tech. (Hangzhou) Co., Ltd.




Building 4, No. 1418-50, Moganshan Road, Gongshu District,
Hangzhou, 310011 Zhejiang, P.R. China



contact@diareagent.com
www.diareagent.com
400-853-5577

Rotavirus and Adenovirus Combo Rapid Test is a rapid visual immunoassay for the qualitative presumptive detection of rotavirus and adenovirus in human fecal specimens.



Simple operation

No special equipment required

Fast results

Clear answer in 10 minutes

High quality combo test

One sample, Two test results

Room temperature stable

2-30°C valid for 24 months

Product Information

Product	REF	Specimen	Test Time	Test Format	Kit Size	Qualification
Rotavirus and Adenovirus Combo Rapid Test	MI-F42001	Human fecal	10 minutes	Multi-Window Device	25 pcs/box	CE

3 Easy Test Steps


1 Prepare Specimen collection and pre-treatment



2 Test Add 2 drops of solution into specimen well




3 Read



Positive



Negative



Invalid

Assure Tech. (Hangzhou) Co., Ltd.

Building 4, No. 1418-50, Moganshan Road, Gongshu District, Hangzhou, 310011 Zhejiang, P.R. China



contact@diareagent.com
www.diareagent.com
400-853-5577